MD

LUPA LED H.F

MIMSAL

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LED lamp with circular magnifying glass.

Head protector with hinge mechanism to protect against dust and scratching.

1. Installation

- Requires a MIMSAL mount, which is installed before the lamp is placed in position.
- A table clamp mount is supplied and wall or ceiling mounts or a rotating base can be supplied as optional accessories.
- MIMSAL cannot guarantee the life and functionality of the product where third party accessories are used.

2. Caution

- Always disconnect the lamp from the power supply when cleaning or carrying out maintenance work.
- The magnifier lamp can cause fire by direct solar radiation. Always pay attention and cover the lens, particularly in areas with direct solar radiation.
- Clean either dry or using a weak solution of normal soap and water. Avoid using alcohol, acids or any type of organic solvent.
- Do not leave on for longer than 8 hours at a time.

3. Troubleshooting and maintenance

- **Maintenance:** Only for qualified electrical engineers or other qualified personnel.
- **Changing the circular LED module:**
 - Always disconnect the lamp from the power.
 - Loosen screws (1 and 3).
 - Remove LED protector (2)./ Change the lens.

Notify by the user / patient to MIMSAL and the competent authority of the country, any serious incident that occurs in relation to the use of the equipment MIMSAL only accepts responsibility for maintenance, repair or changes carried out by MIMSAL staff using original MIMSAL parts.

The manufacturer reserves the right to change the characteristics of its products without prior notice.

	LUPA LED H.F
Light bulb	60 LED H.F 6W
Diopters	5 Dp.
Diameter	12,7 cm
Electrical Protection Class	Class II
International Protection	20
Useful life	40.000 h
Lux at 30 cm	2.000 Lux
Colour temperature	6.500° K
Colour Rendering Index	80%
Length of Arm	90 cm
Color	White
Weight	3 Kg
Power	110/230 V
Frequency	50/60 Hz

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