

INDEX	
Diamond microdermabrasion and sessions	2
Equipment, applicators and full microdermabrasion treatment with F-332	3
Direction of application and diamond tips	4
Diamond microdermabrasion Weelko devices	5
Multifunctional beauty columns and technical sheet	6
Web information	7

## Diamond Microdermabrasion with **MultiEquipment** from Weelko

Acne? **Sun damage?** Blotchy skin? **Fine scars?** Slight wrinkles? **Keratose?** Large pores? **Non-inflamed whiteheads?** Uneven skin texture and glow? **Age spots?** Oil-congested skin? **Skin pigmentation?**

### Light up your skin with **MultiEquipment**, F-332

Modern instrument that combines 4 functions in a single machine: diamond microdermabrasion for a non-invasive and controlled exfoliation treatment; hot and cold hammer to smooth skin after microdermabrasion and improve skin product absorption; facial exfoliation spatula to remove remaining dead cells and ultrasound device to improve cellular metabolism and provide a final healing effect.

Special emphasis is put on microdermabrasion function. Discover everything here inside.



## The diamond microdermabrasion

### About diamond microdermabrasion

It is a non-invasive mechanical exfoliation technique that reduces fine lines and wrinkles by removing outer layers of dead skin. It improves the texture, tone and lightness of the skin and restores its elasticity.

### What does it consist of?

The unit is equipped with an applicator where several diamond tips are adapted. The beautician proceeds with slow sliding motions along the skin surface as shown in page 4. The dead cells are vacuumed back into a waste filter. The heads have different sizes and coarseness for different skin types, areas and depths of polishing.

### Which are the effects?

- Helps natural cell regeneration
- Boosts collagen and elastine production
- Improves skin elasticity
- Reduces the appearance of wrinkles
- Remove acne scars
- Clears up blackheads
- Minimises the size of the pores
- Reduces sun damage

### What is it appropriate for?

- Acne
- Sun damage
- Blotchy skin
- Fine scars
- Light wrinkles
- Keratose
- Large pores
- Non-inflamed whiteheads
- Uneven skin texture and glow
- Age spots
- Congested oily skin
- Skin pigmentation

### Which are the contraindications?

- Inflamed areas
- Bacterial infections
- Cuts and wounds
- Do not suck on eyelids
- Skin tumours
- Freckles
- Angioma
- Cold-sore
- Warts
- Do not use on lips

## Sessions

The recommended time for a simple microdermabrasion session takes around 20-30 minutes. The treatment can be repeated from 3 to 10 times. Treatments should be carried out once a week or once every 15 days. Results are visible from the very first session. It is advisable to perform maintenance sessions every month or every two months.

Treatment	Frequency (days)	Total of sessions
Facial exfoliation	6-10	3-5
Fine lines	7-13	6-8
Neck	12-14	3-5
Acne scars	7-12	8-10
Hyperpigmentation	7-12	8-10

### F-332: Equipment and applicators



**Diamond microdermabrasion :** Exfoliation device with 9 diamond heads and 2 applicators.

**Cold and Hot hammer:** The hot side of the hammer is used to relax muscles and to improve blood circulation. It opens the pores to facilitate the product penetration into the skin. The cold end reduces wrinkles, calms down the skin after microdermabrasion treatments and is used to close pores. Adjustable temperature from 5 to 40°.

**Exfoliation spatula:** Used to exfoliate the outer dermal layer, to treat acne and to remove dead cells, dirt and everyday pollution of the skin. There are 3 types of waves (continuous, pulsed and combined). An exfoliation product must be applied.

**Ultrasound probes:** The ultrasound function is used to stimulate skin tissue and to lift it. This results in a healthier look. It also generates heat below the skin to increase blood flow, dilate capillaries and relax muscles.

#### Full microdermabrasion treatment using F-332:

In order to carry out a complete microdermabrasion treatment with F-332, follow the steps below

**Step 1:** Steam the skin to soften the horny layer and apply a cleansing gel. Then, let the skin dry.

**Step 2:** Carry out the diamond microdermabrasion treatment.

**Step 3:** Apply the cold end of the hammer to calm down the skin surface after the microdermabrasion treatment. Even though it is a non-invasive technique, patient may feel a bit of discomfort after the treatment.

**Step 4:** Use the exfoliation spatula to remove the dead skin that may remain on skin after microdermabrasion treatment.

**Step 5:** Cleanse the skin and apply some gel for ultrasound application.

**Step 6:** Apply the ultrasound probes on the face for a healing effect and to finish the treatment.

**Step 7:** Apply moisturiser with sun protection after the treatment.

**Direction of application:**




The direction of the treatment should be as indicated in the picture, from the central part of the face to the side, from down to up.

- From chin to ears
- From chin to nose
- From nose to temple
- From eye corner to temple
- From central forehead to hairline.

**Diamond tips:**

**G120 and G180:** Ø 22 mm. Body tip, specially for acne, knees, elbows and marks.

**G120, G180 and G240:** Ø 15 mm. Scars, fine lines and skin imperfection. Soft polishing.



**Diamond tips:**

**G180 and G240:** Ø 19 mm. Face. Big spots, large pores, premature ageing, sun exposure.

**G200 and G240:** Ø 8 and 6 mm. Eyes. Both sides of nose. Wrinkles, blackheads, acne, etc.



**Diamond microdermabrasion Weelko devices:**

**U-TECH (F-834):**



**B-EQUIPMENT (F-319A):**



**MULTIEQUIPMENT (F-332)**



**MULTIEQUIPMENT (F-334)**



**MULTIEQUIPMENT (F-336)**



**MULTIEQUIPMENT (F-336B)**



### Multifunctional beauty columns:

F-332 can be used as a standalone product to be placed on any beauty trolley or table. However, it can also be fit on Weelko's beauty columns. **KPSULE** is a tower that allows the storage of 2 devices from the **MultiEquipment** series and it is also equipped with a steamer, a magnifying lamp and a complete trolley to store beauty devices and electrodes. It can also be fit on **MULTILAB**, a trolley with 4 drawers and 2 side doors to store beauty accessories.



### Technical Specifications

Dimensions and weight	Packing size (cm)	49 x 45 x 19,5 cm
	Product size (cm)	33,5 x 9,5 x 22,5
	Gross weight (kg)	5
	Net weight (kg)	4,5
Working parameters		To 40 Hz
Technical specifications	Voltage	220V—240V
	Frequency	50 Hz/60Hz


### Spare parts

No.	PRODUCT	REFERENCE
1	Ultrasound Body probe	10103
2	Ultrasound Face probe	10104
3	Ultrasound Eye probe	10105
4	Diamond head set	40010
5	Foam and metallic filters	40010-1
6	Dermabrasion handle	40016
7	Skin scrubber probe	10018
8	Hot & cold hammer	60050




Visit our TECHNICAL INFO and download all the technical information that you may need. Pictures, dossiers, technical sheets, etc. Do not miss them.

Dossiers | Technical Zoom



Turn on the WEELKO TV in weelko.com. The best way to learn how to use our equipment. Because a picture is worth a hundred words!



weelko.com

Barcelona  
Pol. Ind. La Valldan  
Camí de Garreta,92  
08600 Berga  
Barcelona – Spain  
T. +34 93 824 90 20  
E. [info@weelko.com](mailto:info@weelko.com)