

LS LED

LS LED is an examination lamp. The luminaire is mounted with a highly precise articulated arm and it is manufactured with a metal structure. It is specially designed to be used in medical examination procedures.

Easy to use

Adjustable height articulated arm, allows easy adjustment in different positions, provides a long reach and maximum vertical and horizontal movement, and always remains in the correct position.

Switch located on the connection cable.

Simple and durable

Metallic structure, manufactured with durable materials and with very low maintenance costs.

Mounting Options

The standard version includes AH table clamp.

As for accessories, we offer a wide selection of supports such as wall supports, wall extension arms, trolleys stand and rail brackets.



LS LED

Technical Details	(Cod. LS010316LED)
Light source	LED 7,5W
Head diameter	21 cm
Socket	E-27 ceramic
Arm length	90 cm
Colour	White



Wall support
Cod. B0101800




AH table clamp
Cod. B010190



Rail Plus bracket
Cod. 82026



Extension wall arm
Cod. B010150



Trolley stand
4 kg Cod. 09655