

# LUXIFLEX SENSOR

Luxiflex Sensor is a cold light luminaire, versatile and with a compact design. Luxiflex Sensor is manufactured with a contactless switch. This luminaire is designed for examination procedures that require great performance.

## Flexibility

Flexible arm allows to adjust the height to assure the proper working distance.

## Compact design

Luxiflex offers easy use, thanks to its compact LED head with a diameter of 8 cm and its contactless switch located at the base of the lamp. These features make this model an ideal luminaire for the most demanding examinations procedures.

## Mounting options

The standard version includes AH table clamp.

As for accessories, we offer a wide selection of supports such as wall supports, wall extension arms, trolleys stand and rail brackets.

## High performance LED

Luxiflex Sensor 15.000 lux at 50 centimeters.

## Contactless switch

Allows to switch the lamp on and off without contact.



# LUXIFLEX SENSOR

## Technical Details (Cod. I0104AOLEDSR)

Light source	LED 6W 15°
Lamp lifetime	50.000 hours
Lux a 50 centimeters	15.000 lux
Colour temperature	4.500 °K
Colour rendering	80%
Total length	85 cm
Colour	White



**Trolley stand**  
 4 kg Cod. 09655  
 8,8 kg Cod. 09651



..... Contactless switch



**Wall support**  
 Cod. B0101800



**AH table clamp**  
 Cod. B010190



**Rail Plus bracket**  
 Cod. 82026



**Extension wall arm**  
 Cod. B010150