## MICROSCOPES & REFRACTOMETERS 2023 MICROSCOPES



Digital WiFi microscope KERN ODC-9







Innovative hand-held microscope for mobile applications with immediate display of the image on a smartphone or tablet

## Features

- The digital WiFi hand-held microscope is designed for rapid and simple surface observations. Ideally suited for coins, bank notes, stamps, circuit boards, plants, insects, gems and skin samples for industrial use, for all hobby scientists, children and students
- The KERN ODC 910 WiFi microscope has been specially developed for direct connection to your WiFi-enabled smartphone or tablet with iOS or android
- During the live transfer to your smartphone or tablet you can take photos and videos of the sample you are investigating, and these can also be stored on your device. For larger videos you can also insert a mini SD card directly into the microscope
- With the WiFi microscope you can easily adjust the magnification to suit all conventional samples. The focus can be adjusted to a magnification of 10× as well as 200×

- The six LEDs fitted in a ring shape ensure strong and effective illumination of your sample. Use the adjustment wheel on the microscope to control the illumination setting.
- You can download the app for the ODC 910 WiFi microscope from the Apple App Store or the Android Google Play Store free of charge and this app enables you to directly transfer images and videos from the microscope to your smartphone or tablet through a simple connection
- The scope of delivery includes the WiFi microscope with integrated rechargeable battery pack, a flexible column which is easy to adjust and which has a swan neck so that you can achieve the ideal height setting, as well as a mains adapter

STANDARD WIFI SD RECHARGE 230 V 1 DAY									
Model	Resolution	Interface	FPS	Sensor	Sensor size	Supported operating system	Magnifica- tion levels	Focusing stand	Illumination
KERN									
ODC 910	2 MP	WiFi, SD	15 – 30	CMOS	1/4"	Android, iOS	10×, 200×	Goose neck	6× LED

## **MICROSCOPES & REFRACTOMETERS 2023**

**KERN PICTOGRAMS** 



360° rotatable microscope head



Monocular Microscope For the inspection with one eye



**Binocular Microscope** For the inspection with both eyes



Trinocular Microscope For the inspection with both eyes and the additional option for the connection of a camera



Abbe Condenser With high numerical aperture for the concentration and the focusing of light



Ð

LED

Halogen illumination For pictures bright and rich in contrast

LED illumination Cold, energy-saving and especially long-life illumination



Incident illumination For non-transparent objects



**Transmitting illumination** For transparent objects



Fluorescence illumination For stereomicroscopes

Fluorescence illumination for compound microscopes With 100 W mercury lamp and filter



Fluorescence illumination for compound microscopes With 3W LED illumination and filter



Phase contrast unit For a higher contrast



Darkfield condenser/unit For a higher contrast due to indirect illumination



To polarise the light











 $\infty$ 

INFINITY

Zoom magnification For stereomicroscopes

Infinity corrected optical system

Infinity system

Auto-focus





Parallel optical system For stereomicroscopes, enables fatigue-proof working

For automatic control of the focus level



SD card For data storage



USB 2.0 digital camera For direct transmitting of the picture to a PC USB 2.0



USB 3.0 digital camera For direct transmitting of the picture to a PC



WIFI data interface: For transmitting of the picture to a mobile display device



HDMI digital camera For direct transmitting of the picture to a display device



PC software To transfer the measurements from the device to a PC.



Automatic temperature compesation For measurements between 10  $^\circ\mathrm{C}$  and 30  $^\circ\mathrm{C}$ 



Protection against dust and water splashes IPxx:

The type of protection is shown in the pictogram cf. DIN EN 60529:2000-09, IEC 60529:1989+A1:1999+A2:2013



C-Mount	Adapter for the connection of a camera to a trinocular microscope			
FPS	Frames per second			
H(S)WF	High (Super) Wide Field (Eyepiece with high eye point for wearers of glasses)			
LWD	Long Working Distance			
N.A.	Numerical Aperture			
SLR camera	Single-Lens Reflex camera			
SWF	Super Wide Field (Field number at least Ø 23 mm for 10× eyepiece)			
W.D.	Working Distance			
WF	Wide Field (Field number up to Ø22 mm for 10× eyepiece)			
	······ ·······························			



operation

Battery operation Ready for battery operation. The battery type is specified for each device.



Battery operation rechargeable Prepared for a rechargeable battery

230 V

Plug-in power supply 230V/50Hz in standard version for EU. On request GB, AUS or USA version.



Integrated power supply unit Integrated in microscope. 230V/50Hz standard EU. More standards e.g. GB, AUS or USA on request.



Package shipment The time required to manufacture the product internally is shown in days in the pictogram.