

MAG2000

Mag2000 is the world's best-selling low-frequency PEMF device.

It is ideal for the treatment of various bone and muscle diseases thanks to its 35 preset programs.



TECHNICAL FEATURES

- ✓ 35 programs
- ✓ 2 output channels
- ✓ Intensity up to 150 Gauss for single channel
- ✓ Frequency from 1 to 100 Hz

MEDICAL DEVICE CLASS IIa
MAG2000 is a medical device CE0068.

KIT

- ✓ MAG2000
- ✓ Elastic therapeutic belt with solenoids
- ✓ Medical power adaptor
- ✓ User Manual
- ✓ 1 PVC envelope
- ✓ 1 magnet tester
- ✓ Bag



DIMENSIONS 179 x 107 x 50 mm

MAG2000

MAG2000 has 20 pre-setted programs for the treatment of several pathologies, 14 programs identified by the working frequency (from 1 to 100 Hertz) and an Autoscan program with continuous frequency variation for the treatment of both hard and soft tissues at the same time.

MAG2000 has the possibility to set the intensity of magnetic field for all the programs up to 300 Gauss on 2 channels (150 Gauss for single channel).

PROGRAMS

01 Osteoporosis	17 Coxartrosis	33 90 Hz
02 Arthrosis	18 Muscular Atrophy	34 100 Hz
03 Arthritis	19 Muscular Contracture	35 Autoscan
04 Cervical arthrosis	20 Osteonecrosis	
05 Articular pain	21 1 Hz	
06 Cervicalgias	22 3 Hz	
07 Sprains	23 5 Hz	
08 Fractures	24 10 Hz	
09 Epicondylitis	25 15 Hz	
10 Epitrocleitis	26 20 Hz	
11 Intercostal contusions	27 30 Hz	
12 Lumbalgy	28 40 Hz	
13 Lumbar pain	29 50 Hz	
14 Shoulder arthrosis	30 60 Hz	
15 Knee arthrosis	31 70 Hz	
16 Periarthritis	32 80 Hz	