


MEASURING MAT FOR BABIES

seca

Although rugged, the surface is easy to care for and can be cleaned with all commercially available disinfectants




Size 125x14x30 cm

27329

• **27329 SECA 210 - BABY MEASURING MAT**
Mobile measuring mat for babies and toddlers. Kind to the skin, washable, foldable. The space saving solution for easy and precise measurement of baby length. Range 10-99 cm. Graduation: 5 mm.

CE 0123



CE 0476

• **27331 GIMA BABY MEASURING MAT**
Washable and foldable measuring mat for babies and toddlers, easy to carry. For easy and precise measurement of baby length. Size: 121x30 cm. Measuring range: 10-99 cm. Graduation: 5 mm.

SECA INFANTOMETER

CE 0123

The unique V-shape automatically places the baby in the right measuring position

seca



The head and foot positioners lock into place to capture results. Indication can be read after measuring

• **27316 SECA 416 - INFANTOMETER**

Infantometer for measuring babies and toddlers in hospital. The sturdy construction and stability of the infantometer makes it quick and easy to measure babies and toddlers up to two years old. The board is generously designed and the raised and softly rounded sides guide the baby easily and safely into the right position. Permanently mounted with rollers on two guide rails, the footpiece slides smoothly along the scale. Because the footpiece is lockable, the baby can be removed and the result noted later. All parts are tough and durable. Measuring range: 33-100 cm. Graduation: 1 mm. Size: 110x16.9x40.2 cm. Weight: 3.8 kg.

GIMA PEDIATRIC MEASURING INSTRUMENTS

CE 0476


27326



- **27326 BABY MEASURE "CALIBRO"**
Made of oxidized aluminium, metric graduation 0-900 mm, division 5 mm. Certified by Italian "weights and measure institute". Weight: 300 g.
- **27327 CRANIOMETRO - HEAD MEASURE**
Made of special, soft, non stretchable plastic. Measures child head girth up to 60 cm.
- **27337 ORCHIDOMETER - wood**
Used by pediatricians for screening of children's testicles

CE 0476

27327



27337

A string of 12 numbered wooden beads of increasing size from about 1 to 25 mm used to determine size of testicles.

HEIGHT MEASURE

- **27335 TAPE HEIGHT MEASURE**
Clear scale 0-200 cm. Easy wall fixation. Supplied with screw for wall mounting.


CE 0476



ELECTRONIC HEIGHT ROD

SOEHNLE PROFESSIONAL
LCD digit height: 9 mm

27333



Size: 345x92x28 mm Power supply: 2x1.5 V AA batteries Weight: 300 g

- **27333 SOEHNLE ELECTRONIC HEIGHT ROD**
Measuring range: 50-240 cm (20-94.5 inches), division: 1 cm. Manual in: DE, GB, FR, ES, IT, GR, CZ.

Switchable to inches

Size: 630x290x70 mm

Power supply: 2x1.5V AAA batteries

- **27264 SOEHNLE BABY HEIGHT ROD**
Digital and mechanical rod 35 to 80 cm, division 1 cm. LCD display, digit height 16 mm. Manual in 23 languages.

27264



BONE MEASURE

- **27325 BONE MEASURE**
Slide-knot caliper in anticorodal to measure knee, shoulder, wrist, elbow... Measuring capacity: 150 mm. Burner length: 100 mm. Tested by Italian "weights and measures institute".

CE 0476

27325

