



GIMA



CE 0476

Via Marconi, 1 - 20060 Gessate (MI) Italia
Web Site: www.gimaitaly.com

IP 21

SUPERVEGA EVO

REF 28183

LOT XXXXXXXXX

SN XXXX



14 V \approx 4A tramite alimentatore esterno o alimentazione interna (batteria Li-ion 14,8V \approx 5,2A) o con adattatore accendisigari 12V \approx 4A
14V \approx 4A by external power adaptor or internal powered equipment (Li-ion battery 14,8V \approx 5.2A) or with cigarette lighter adaptor 12V \approx 4A



Gima S.p.A. - Via Marconi, 1
20060 - Gessate (MI) - Italy
Made in Italy



Quando alimentato tramite batteria dell'auto o adattatore AC/DC - When powered by car battery adaptor or AC/DC adaptor



Rev.1-01.21