



1pz.



REF 34337

LOT 0000xx

SET LARINGO "MAXLITE" F.O. LED 4 LAME MC-INT 1-2-  
MAXLITE F.O.LED ADULT SET 4 BLADES MC-INTOSH 1-2-3  
SET LED À F.O. MAXLITE ADULTE 4 lames: Mc Intosh 1-2-3  
SET LARINGO "MAXLITE" F.O. LED ADULTOS 4 HOJAS MC-INT 1-2-3



GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy  
Made in Pakistan - [www.gimaitaly.com](http://www.gimaitaly.com)



509037/C - Rev. 2 - 04/20  
Size 19x23,5x4,8 cm

- REF 34333** Adult set 3 blades Maxlite MCInosh 1-2-3 / Xenon/Helogen bulb (2.5V)  
**REF 34334** Adult set 4 blades Maxlite MCInosh 1-2-3-4 / Xenon/Helogen bulb (2.5V)  
**REF 34335** Adult set 3 blades Maxlite 2-3-4 / Led bulb (2.5V)  
**REF 34337** Adult set 4 blades Maxlite MCInosh 1-2-3-4 / Led bulb (3.5V)  
**REF 34338** Adult set 4 blades Maxlite MCInosh 1-2-3-4 / Lampadina Led (3.5V)
- REF 34332** Pediatric set 3 blades Maxlite Miller 00-0-1 Xenon/Helogen bulb (2.5V)  
**REF 34336** Pediatric set 3 blades Maxlite Miller 00-0-1 Led bulb (2.5V)
- Lampadina Xenon/Allogena (2.5V)  
 Set pediatrico 3 lame Maxlite Miller 00-0-1  
 Lampadina Led (2.5V)

## Laringoscopi a Fibre Ottiche Maxlite




- Laryngoscope à fibre optique - Maxlite
- Laringoscópios em fibra ótica - Maxlite
- Laryngoskopy światłowodowe - Maxlite
- Λαρυγγοσκόπιο οπτικής ίνας - Maxlite



- Luce fredda ad elevata luminosità - fibra Ø 4,5 mm
- Autoclavabile a 134°C per 2000 cicli
- Intensità luminosa 8000 LUX (2,5V) o 14000 LUX (3,5V)
- Fascio di fibre ottiche integrato, senza cavità che potrebbero trattenere sporcizia o fluidi corporali.

[www.gimaitaly.com](http://www.gimaitaly.com)



Gima S.p.A.  
Via Marconi, 1 - 20060 Gessate (MI) Italy  
[gima@gimaitaly.com](mailto:gima@gimaitaly.com) - [export@gimaitaly.com](mailto:export@gimaitaly.com)  
[www.gimaitaly.com](http://www.gimaitaly.com)  
Made in Pakistan

Vedere etichetta bar code  
See bar code label



**GIMA**  
Via Marconi, 1 - 20060 Gessate (MI) Italy  
[gima@gimaitaly.com](mailto:gima@gimaitaly.com) - [export@gimaitaly.com](mailto:export@gimaitaly.com)  
[www.gimaitaly.com](http://www.gimaitaly.com)

## Maxlite Fiber Optic Laryngoscopes



- Laringoscopios de fibra óptica - Maxlite
- Maxlite LWL-Laryngoskope
- Laringoscop cu fibră optică - Maxlite
- مناظير الحنجرة بألياف بصرية - Maxlite



- Superior cold, bright light Ø 4,5 mm fibre
- Autoclavable at 134°C for 2000 cycles
- Light intensity 8000 LUX (2.5V) or 14000 LUX (3.5V)
- Integrated F.O. bundle with no cavities to trap dirt or body fluids

[www.gimaitaly.com](http://www.gimaitaly.com)



**GIMA**  
Via Marconi, 1 - 20060 Gessate (MI) Italy  
[gima@gimaitaly.com](mailto:gima@gimaitaly.com) - [export@gimaitaly.com](mailto:export@gimaitaly.com)  
[www.gimaitaly.com](http://www.gimaitaly.com)