



CE 0476

Via Marconi, 1 - 20060 Gessate (MI) Italia
Web Site: www.gimaitaly.com

IP 22

SUPERVEGA EVO

REF 28182

LOT XXXXXXXXX

SN XXXX



Alimentazione interna (Batteria li-Ion
14,8V 5.2A) o con adattatore per
ambulanza (12V 4A)

Internally powered equipment (li-Ion Battery
14,8V 5.2A) or with Ambulance
adapter (12V 4A)

Gima S.p.A. - Via Marconi, 1
20060 - Gessate (MI) - Italy
Made in Italy

Quando alimentato dall'adattatore per
ambulanza / When powered by
Ambulance adapter



Rev.1-01.21