



CLINNOVA Clinical 2 Section



A very strong and rigid couch available with a wide range of accessories. Suitable for supporting very heavy patients for examination and minor treatment procedures.

Models

See Page 2 for all available models



Model NV2550-CLS
Colour dark grey



Model NV2560-PRM-EPFS
Shown with optional perimeter foot switch
Colour Lotus Green



Negative head section to -25°

Standard Features

- Safe working load and lifting capacity of 265Kg (583lbs)
- Gas-assisted backrest adjustment from -20° to +80° on Models NV2550 & NV2560 with anti-entrapment design
- Electric backrest adjustment on NV2570
- Paper Roll Holder fitted as standard
- Choice of hydraulic or electric height elevation
- Electric models raise from low to high in only 15 seconds
- Excellent height range from 49cm to 99cm
- Linked retractable wheel design with adjustable foot
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Anti-bacterial powder coated frame
- 70cm wide Easy Clean Upholstery as standard
- Choice of 15 standard upholstery colours

Optional Accessories

- 6004** Breathing hole & plug
- 6200** Single foot switch for NV2560
- 6207** Double foot switch for NV2570
- 9002** Side support rails (Pair)
- 6002P** Side support rail cushions (Pair)
- 6007** Removable head/knee support
- 6057** Battery back-up facility

Models Available



NV2550-CLS

Hydraulic Elevation
Gas Assisted Backrest
Classic Base (CLS)



NV2560-CLS

Electric Elevation
Gas Assisted Backrest
Classic Base (CLS)



NV2570-CLS

Electric Elevation
Electric Backrest
Classic Base (CLS)



NV2550-PRM

Hydraulic Elevation
Gas Assisted Backrest
Premium Base (PRM)



NV2560-PRM

Electric Elevation
Gas Assisted Backrest
Premium Base (PRM)



NV2570-PRM

Electric Elevation
Electric Backrest
Premium Base (PRM)