


## CLINNOVA Mobile Hygiene Tables



For environments such as special schools, where many of the children require regular hygiene attention, these variable height tables provide not only a low transfer height but quick and simple adjustment to the correct height for changing procedures, reducing back and postural problems for support staff. Medium & Large Tables are also suitable for adults & juniors that require additional care. These tables are fitted with 125mm castors and side support rails as standard.

### Models

See Pages 2 & 3 for all available models



### Standard Features

- Safe working load and lifting capacity of 265Kg (583lbs)
- Choice of hydraulic or electric height elevation
- Choice of Small, Medium or Large upholstery lengths
- Choice of individually braking castors (IBC) or large mobile wheel design (LMWD) with central locking castors and steering facility
- Excellent variable height range from 46 to 99cm
- Electric models raise from low to high in only 15 seconds
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Anti-bacterial powder coated frame
- 70cm wide Easy Clean Upholstery as standard
- Choice of 15 standard upholstery colours

### Optional Accessories

- 6200** Single Foot Switch (*Electric only*)
- 6002P** Side Support Rail Cushions
- 6007** Head/Knee Support Cushion



Optional 6002P Cushions for Side Support Rails



Model ST1651A-PRM (Large)  
Colour Olive



Model ST1461B-CLS (Medium)  
Colour Burgundy



Model ST1661C-PRM (Small)  
Colour Sky Blue

## Models Available



### NV1451A-CLS

Hydraulic Elevation  
Large  
Independently Braking Castors  
Classic Base (CLS)



### NV1461A-CLS

Electric Elevation  
Large  
Independently Braking Castors  
Classic Base (CLS)



### NV1451B-CLS

Hydraulic Elevation  
Medium  
Independently Braking Castors  
Classic Base (CLS)



### NV1461B-CLS

Electric Elevation  
Medium  
Independently Braking Castors  
Classic Base (CLS)



### NV1451C-CLS

Hydraulic Elevation  
Small  
Independently Braking Castors  
Classic Base (CLS)



### NV1461C-CLS

Electric Elevation  
Small  
Independently Braking Castors  
Classic Base (CLS)



### NV1451A-PRM

Hydraulic Elevation  
Large  
Independently Braking Castors  
Premium Base (PRM)



### NV1461A-PRM

Electric Elevation  
Large  
Independently Braking Castors  
Premium Base (PRM)



### NV1451B-PRM

Hydraulic Elevation  
Medium  
Independently Braking Castors  
Premium Base (PRM)



### NV1461B-PRM

Electric Elevation  
Medium  
Independently Braking Castors  
Premium Base (PRM)



### NV1451C-PRM

Hydraulic Elevation  
Small  
Independently Braking Castors  
Premium Base (PRM)



### NV1461C-PRM

Electric Elevation  
Small  
Independently Braking Castors  
Premium Base (PRM)

## Models Available



### NV1651A-CLS

Hydraulic Elevation  
Large  
Central Locking Castors  
Classic Base (CLS)



### NV1661A-CLS

Electric Elevation  
Large  
Central Locking Castors  
Classic Base (CLS)



### NV1651B-CLS

Hydraulic Elevation  
Medium  
Central Locking Castors  
Classic Base (CLS)



### NV1661B-CLS

Electric Elevation  
Medium  
Central Locking Castors  
Classic Base (CLS)



### NV1651C-CLS

Hydraulic Elevation  
Small  
Central Locking Castors  
Classic Base (CLS)



### NV1661C-CLS

Electric Elevation  
Small  
Central Locking Castors  
Classic Base (CLS)



### NV1651A-PRM

Hydraulic Elevation  
Large  
Central Locking Castors  
Premium Base (PRM)



### NV1661A-PRM

Electric Elevation  
Large  
Central Locking Castors  
Premium Base (PRM)



### NV1651B-PRM

Hydraulic Elevation  
Medium  
Central Locking Castors  
Premium Base (PRM)



### NV1661B-PRM

Electric Elevation  
Medium  
Central Locking Castors  
Premium Base (PRM)



### NV1651C-PRM

Hydraulic Elevation  
Small  
Central Locking Castors  
Premium Base (PRM)



### NV1661C-PRM

Electric Elevation  
Small  
Central Locking Castors  
Premium Base (PRM)