



Width  
70cm



## CLINNOVA Therapy 3 Section – Standard



A very strong and rigid 3 section couch with a traditionally long head section to suit a wide range of physiotherapy and general examination procedures. The couch can be adjusted to provide support for the patient in prone, supine and seated positions. The couch has an excellent height range of 46cm to 95cm and is fitted with a breathing hole & plug as standard

### Models

See Page 2 for all available models

Model NV3567-CLS  
Colour Red



### Standard Features

- Safe working load and lifting capacity of 265Kg (583lbs)
- Choice of hydraulic or electric height elevation
- Gas-assisted backrest adjustment from 0° to +85°
- Gas-assisted footrest adjustment from -20° to +80°
- Electric models raise from low to high in only 15 seconds
- Excellent height range from 46cm to 95cm
- Linked retractable wheel design with adjustable foot
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Anti-bacterial powder coated frame
- 70cm wide Easy Clean Upholstery as standard
- Choice of 15 standard upholstery colours

### Optional Accessories

- 6200** Single Foot Switch
- 9003** Side Support Rails
- 6003P** Side Support Rail Cushions
- 9005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 9001** Electric Perimeter Foot Switch (EPFS)



Offset gas strut for unobstructed patient view



Negative head/foot section to -25°  
Shown with optional paper roll holder



Low patient transfer height  
Shown with optional paper roll holder

## Models Available

---



### **NV3557-CLS**

Hydraulic Elevation

Classic Base (CLS)



### **NV3567-CLS**

Electric Elevation

Hand or Foot Switch

Classic Base (CLS)



### **NV3567-CLS-EPFS**

Electric Elevation

Perimeter Foot Switch

Classic Base (CLS)



### **NV3557-PRM**

Hydraulic Elevation

Premium Base (PRM)



### **NV3567-PRM**

Electric Elevation

Hand or Foot Switch

Premium Base (PRM)



### **NV3567-PRM-EPFS**

Electric Elevation

Perimeter Foot Switch

Premium Base (PRM)

## Models Available

---



**NV3557-CLS**

Hydraulic Elevation  
Classic Base (CLS)



**NV3567-CLS**

Electric Elevation  
Hand or Foot Switch  
Classic Base (CLS)



**NV3567-CLS-EPFS**

Electric Elevation  
Perimeter Foot Switch  
Classic Base (CLS)



**NV3557-PRM**

Hydraulic Elevation  
Premium Base (PRM)



**NV3567-PRM**

Electric Elevation  
Hand or Foot Switch  
Premium Base (PRM)



**NV3567-PRM-EPFS**

Electric Elevation  
Perimeter Foot Switch  
Premium Base (PRM)