



Innovation Deluxe 2 Section



The Innovation Deluxe 2 section couch is suitable for use within specialist critical care environments, endoscopy, ultrasound scanning and minor treatment procedures. Electric height, backrest and 2-way tilt ensure smooth transition when positioning patients for specialist procedures and treatments. Auto-electric CPR functionality and auto-levelling software return the couch quickly to a fully horizontal resuscitation position in the event of an emergency.

Models

- MG2675** Electric with electric backrest and 2-way tilt. Large central locking wheel design with steering facility
- MG2475** Electric with electric backrest and 2-way tilt. Large independently locking wheel design



Model MG2675
with large central locking wheel
design with steering facility

Deluxe Features

- Electric height range from 47cm to 97cm
- Electric trendelenburg tilt to $-/+15^\circ$
- Electric profiling backrest angle from -20° to $+85^\circ$
- Safe working load and lifting capacity of 250Kg (550lbs)
- Auto-electric CPR function via handset
- Auto-intelligent operation
- Vertical lifting design to minimise space and provide maximum stability throughout all profiling positions
- Smooth moulded base cover minimises dirt traps and facilitates easy cleaning
- Anti-bacterial powder coated frame and base cover
- Large Tente® $\varnothing 125\text{mm}$ castors
- **SEERS ComfortPlus™** Upholstery - features anti-microbial vinyl incorporating pressure relieving design for additional patient comfort with a foam depth of 60mm
- Standard width upholstery of 70cm
- Clearance beneath frame for mobile hoist
- Head support cushion fitted on Deluxe models
- Paper roll holder fitted on Deluxe models
- Choice of upholstery colours, for full range of colour options see page 51

Optional Accessories

- 6003** Side support rails (Pair)
- 6003P** Side support rail cushions (Pair)
- 6004** Breathing hole and plug for prone treatments
- 6006** Upper limb support
- 6054** Telescopic IV pole
- 6057** Battery back-up facility

See page 50-51 for additional information on our range of optional accessories



Negative head section to -20°



Trendelenburg tilt to -15°