



Medicare Bariatric 2 Section



The Bariatric 2 Section Couches feature electric height and backrest operation for effortless patient positioning for various medical examinations and treatment procedures. There is a choice of widths available - Models SM2573 and SM2673 are available with a width of 80cm and Models SM4542 and SM4642 which have a wide frame, are available at 105cm as standard and 125cm at an additional cost. All couches provide a safe working load and lifting capacity of 325Kg (715lbs).

Models

- SM2573** Electric with electric backrest 80cm wide. Retractable wheel system
- SM2673** Electric with electric backrest 80cm wide. Large wheel design, with central locking and steering facilities
- SM4542** Electric with electric backrest 105cm wide. Retractable wheel system
- SM4642** Electric with electric backrest 105cm wide. Large wheel design, with central locking and steering facilities



Model SM2573 - standard width 80cm.
Colour sky blue

Standard Features

- Safe working load and lifting capacity of 325Kg (715lbs)
- Choice of retractable wheel or large wheel design, with central locking & steering facility
- Adjustable angle electric backrest operation from 0° to +75°
- Electric height elevation from 45cm to 98cm
- Purpose made precision bearings on all pivot points
- Anti-bacterial powder coated frame
- Clearance beneath frame for mobile hoist
- Choice of 80cm, 105cm or 125cm width upholstery
- Hand switch fitted as standard
- Choice of upholstery colours, for full range of colour options see page 51

Optional Accessories

- 6076** Memory foam upholstery (60mm) for additional patient comfort
- 6004** Breathing hole and plug for prone treatments
- 6074** Side support rail (Pair)
- 6074P** Side support rail cushions (Pair)
- 6005** Paper roll holder
- 6207** Double foot switch

See page 48-49 for additional information on our range of optional accessories



Model SM4542 shown with optional fold-down side support rails and cushions



Model SM4642 - standard width 105cm with central locking wheels and steering facility. Colour sky blue