

MULTIEQUIPMENT

Ultrasound instrument equipped with 3 probes for face, eyes and body.

Equipo de ultrasonidos equipado con 3 sondas para ojos, cuerpo y cara.

Appareil à ultrasons équipé de 3 sondes pour le visage, les yeux et le corps.

Ultraschallgerät mit 3 verschiedenen Sonden (Orbicular, Gesicht und Körper) ausgestattet.

Apparecchiatura di ultrasuoni dotata di tre sonde (facciale, corporale e orbicolare).

**Ultrasound
Ref. F337**

**Weight: 2,2Kg
Size: 33x27,5x10 cm**

Input voltage: 100~120 /220 ~240V
Input frequency: 50Hz/60Hz
Rated power: 25W
Working frequency: 1MHz

**PACKING SIZE cm
49x45x19,5**

**N.W. kg
2,2**

**G.W. kg
3,3**

