

UV-POWER

Ultra-violet steriliser device conceived for plastic non-sharp instruments.

Esterilizador de luz ultravioleta para objetos de plástico e instrumental no punzante.

Stérilisateur à ultraviolette pour le petit matériel en plastique non pointu.

UV Sterilisor ideal für nicht spitzige Instrumente und auch für Kunststoffgegenstände.

Sterilizzatore a luce ultravioletta per oggetti di plastica e per strumenti non pungenti.

**UV Sterilizer
Ref. S02**

Size: 40x24x21cm

Weight: 2 Kg

Input voltage: 220 ~240V

Input frequency: 50/60Hz

Rated power: 12W

UV lamp lifetime: 5000 hours

Regular working time: > 1 year

PACKING SIZE cm

43x26x22

N.W. kg

2

G.W. kg

3

